## 基因剔除(剔入)細胞株核心實驗室 Gene Knockout/in Cell Line Modeling Core

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**E-CRISPR**  
Design of CRISPR constructs

dkfz. GERMAN CANCER RESEARCH CENTER

De-novo Evaluation Database Forum Help Links

E-CRISPR is a software tool to design and evaluate target sites for use with the Clustered Regularly Interspaced Short Palindromic Repeat (CRISPR) system.

The web application uses fast algorithms to identify sgRNA target sequences in any nucleotide sequence for use in CRISPR/Cas mediated genome editing. E-CRISPR analyzes target specificity of the putative designs and assesses their genomic context (e.g. exons, transcripts, CpG islands). The design process incorporates different options of how CRISPR constructs can be used in experimental applications.

**Recent Updates:**

- 26 May 2014 We are happy to announce a further major update to our E-CRISPR web service. Many new organisms have been added together with big changes in the web front end. Hence you will find the new forum and many other new things here in the new BETA version 3.1.
- 14 April 2014 In this minor update different default values for de-novo sgRNA design have been implemented, allowing for more designs to be found.
- 01 April 2014 The following organisms have been added to E-CRISPR:  
Zea mays  
Ustilago hordei  
Toxoplasma gondii ME49  
Gasterosteus aculeatus  
Populus trichocarpa
- 20 March 2014 A new version of E-CRISPR has been released (version 3.0). It includes more off-target search options and we implemented speed improvements to enable the design of sgRNAs against up to 200 genes in parallel.

**De-NOVO**

```
>seq
ACTGATTATTGGGAGTCGATCGTAGTCGATTAGCTAGCATAATAGCGTACTGAT
TGGCTATGCTGAGCT
```

**EVALUATION**

```
>seq1
ACCGATATCGATGCTGGATATCGATACATCGATGTAGCTAGCTGATC
>seq2
TTAGGGACTGACCGGAGTATAGCGGATTA CGGGCGATTACGGCGATTA
```

The diagram illustrates the workflow from a DNA sequence to the evaluation of target sites. The 'De-NOVO' step shows a sequence with a red arrow pointing to a target site. The 'EVALUATION' step shows the same sequence with a red arrow pointing to a target site and a purple arrow pointing to a target site, indicating the evaluation of different target sites.

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## E-CRISP

Design of CRISPR constructs

dkfz. GERMAN  
CANCER RESEARCH CENTER

De-novo Evaluation Database

Forum Help

Links

### 1. Select organism:

Homo sapiens sapiens (Human, GRCh37.75) [HELP]

### 2. Enter target sequence:

FASTA format sequence(s)

FASTA example [HELP]

```
>TBK
AGCGAGCTCGGCCGCCGCGGAGTACGCACGCGCAGGCGCGCCCCGC
ACTCGGTCACGTGCTCCGTGGCCCGGGCTGGCGAAGCCGGAAGTAG
CCTGGGGCGCGAGAAGGCCCGGGAGCCGCGGGCTGTACGCGCGGA
CACTCGCGGGGTAAAGAGGGGCGGGCGAGGCCGGTGGGGCGGACG
```

Alternatively upload a file in FASTA format or a new line separated list of gene symbols and/or accession numbers

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[HELP]

Clear

### 3. Choose parameter presets:

Relaxed Medium **Strict** Single design

選擇Single design (CRISPR Wt)或  
Paired design (CRISPR D10A)

### 4. Start web application:

Search

Start sgRNA analysis

Display advanced options

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[click here to change CRISPR properties](#)

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bp Minimum guide RNA length after PAM  
[HELP]

bp Maximum guide RNA length after PAM  
[HELP]

% < G <  %

% < A <  %

% < T <  %

% < C <  %

bp 3' flanking sequence length [HELP]

bp 5' flanking sequence length [HELP]

bp Expected flank size (for restriction site analysis)

bp Tagging window downstream of the codon

bp Tagging window upstream of the codon

Number of coding exons downstream the start codon for K.O.

bp Minimal spacer length for paired designs


bp Maximum spacer length for paired designs

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Exclude non-gene hits [HELP]

Exclude non-exon hits [HELP]

Exclude hits outside of coding sequences

Exclude hits in CpG islands [HELP]

Transcript specificity [HELP]

Exon specificity [HELP]

Retrieve and save a recombination donor matrix

Assess restriction sites of the whole sequence [HELP]

Restriction enzyme library [HELP]

#### 7. Off-target analysis:

select bowtie version

select off-target database

Exclude designs with more than X off-targets

Number of 5' mismatch positions

Tolerated edit distance to the target sequence

Bowtie2 pre-sets

Select to check for secondary off-targets [HELP]

Please choose the targets to test and/or paste a custom sequence into the text area below.

The pasted sequence should be FASTA format.

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#### 8. Output:

Maximum number of results per exon

- Produce a GFF formatted output file
- Create an image showing genomic context
- Output the result table to the browser window
- Produce additional information for the Matchstring (Warning: This option might take a while!)

Start

This is a beta version. If you have any questions or find potential bugs please do not hesitate to contact us.  
Go to original E-CRISP

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## E-CRISP

Design of CRISPR constructs

dkfz. GERMAN  
CANCER RESEARCH CENTER

De-novo

Evaluation

E-TALEN

E-RNAi


GenomeRNAi

Help About

While you are waiting, why don't you visit E-RNAi or E-TALEN?

Did you know that:

The **nails** on long fingers grow faster than short fingers; and fingernails grow four times faster than toenails.

This page will be automatically updated every 28 seconds until search is done.

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GERMAN  
CANCER RESEARCH CENTER

De-novo
Evaluation
E-TALEN
E-RNAi
GenomeRNAi
Help
About

[Download a tabular report for all query sequences together](#)

---

**Query name:** TBK **Query length:** 483 **Query location:** ::0::0

Total number of possible designs = 174  
 Number of successful designs = 4  
 Number of designs that hit a specific target = 71  
 Number of designs excluded because they hit multiple targets or none = 37  
 Number of designs excluded because their nucleotide composition was too invariable or contained TTTT = 3  
 Number of designs excluded because the maximum of designs per exon was exceeded = 67

Name	Nucleotide sequence	Ranking	Target	Matchstring	Number of Hits
TBK_6_96	TTCTTCCCCCTTATCG TGAGG	109.52	ENSG00000058668::A TP2B4	<a href="#">Matchstring Info</a>	1
TBK_28_384	AACTAGCACCCTAGG TAGG	100	5	<a href="#">Matchstring Info</a>	2
TBK_28_384	AACTAGCACCCTAGG TAGG	100	ENSG00000145675::PI K3R1	<a href="#">Matchstring Info</a>	2
TBK_21_384	AGTTTCTGAGTTTACA GAATGG	87.27	8	<a href="#">Matchstring Info</a>	3
TBK_21_384	AGTTTCTGAGTTTACA GAATGG	87.27	ENSG00000175497::D PP10	<a href="#">Matchstring Info</a>	3
TBK_21_384	AGTTTCTGAGTTTACA GAATGG	87.27	1	<a href="#">Matchstring Info</a>	3
TBK_1_0	AAACACTATCTCCAGA GGG	70	11	<a href="#">Matchstring Info</a>	5

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


## 方法二、Restriction Enzyme (RE) digestion assay

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